



LP ITALIANA SPA

Via C. Reale 15/4, 20157 Milano, ITALY
T +39 02 3933 06.1 - F +39 02 3931 3484
www.lpitaliana.com - info@lpitaliana.com

Capitale Sociale € 309.600,00
R.E.A. MI 882798
Reg. Imp. C.F. e P.I. MI 01794050151

EC DECLARATION OF CONFORMITY

In accordance with the 98/79/CE – ANNEX III, points applicable to:

Test tubes for samples of human body fluids

References:

- 4th group as per Declaration according to D.Lgs. 332 submitted to the Italian Ministry of Health on 08/05/2001
- DIR 98 / 79 /CE acknowledged by Italy with D.Lgs. 332 of 08/09/2000
- Device Master File Nr. 02-20- "Test tubes for in vitro diagnostic examination"
- Certified Quality Management System to UNI EN ISO 13485:2016

LP ITALIANA SPA

declares

that concerned products from its own catalogue, referred to the product codes in the enclosed list 02-20, fulfil the obligations imposed by the Directive, which apply to them, and

meet

the essential requirements provisions set out in Annex I (which apply to them) and

are not contained

in the List A nor List B of **Annex II** of the Directive.

LP ITALIANA SPA keeps at the disposal of the Competent Authorities this Declaration of Conformity, together with the relative Device Master File (cod. FT 02-20) comprising the following documents.

| | |
|--|------------------------------|
| Technical: <i>General product description including variants</i> | <i>Design information</i> |
| <i>Documentation of the Quality System</i> | <i>Technical drawings</i> |
| <i>Direction for use description</i> | <i>Risk analysis</i> |
| <i>Characteristics of the basic materials</i> | <i>Microbiological state</i> |
| <i>Performances and limits</i> | <i>Tests reports</i> |
| <i>Method of manufacture</i> | <i>Labelling</i> |

and assuring that the relative method of manufacture follows the principles of quality assurance appropriate for these test tubes:

Methods of manufacture according to adequate QA principles
Organizational structure and responsibility
Applicable Operating Procedures of LP ITALIANA SPA's Quality Management System
Applicable Operating Instructions of LP ITALIANA SPA's Quality Management System
Statistical control of production quality
Quality Management System performance control means

LP ITALIANA SPA has instituted and keeps up to date its PO 04.03 systematic procedure, covering the following aspects:

- *Reviewing the experience gained from test tubes in the post-production phase*
- *Implementation of any appropriate means to apply any necessary corrective actions taking account of the nature and risk in relation to the use of the test tubes*
- *Immediate notification to the Competent Authorities of incidents involving the use of the test tubes, on learning of them*

With the above, LP ITALIANA SPA notifies the Italian Ministry of Health to affix on the test tubes, subject of this Declaration, the CE marking of conformity in accordance with article 16 (17) and Annex X of the mentioned Directives.



LP ITALIANA SPA
Legal Representative
(*Francesco Leopardi*)

Signed in Milan the, 13 December 2021



| | | | | |
|-------------------------------|--|------|------------|-----------|
| List Number: | DMF 02-20 Rev. 11 | Date | 2019.02.28 | Sheet 1/4 |
| EC Declaration of Conformity: | "Test tubes for samples of human body fluids" | | | |
| Product classification: | IVD - "Other Device" | | | |
| EC Directive | 98/79/CE | | | |
| Device Master File Nr | FT 02-20 | | | |



Product Code List

| Code | Description  | Descrizione  |
|--------|---|--|
| 116047 | 50 ml Polypropylene test tube, Ø 30 x 115 mm, with blue screw cap, graduated | Provetta in polipropilene da 50 ml, Ø 30 x 115 mm, graduata, con tappo a vite blu |
| 116048 | 50 ml Polypropylene test tube, Ø 30 x 115 mm, with blue screw cap, graduated, sterile R | Provetta in polipropilene da 50 ml, Ø 30 x 115 mm, graduata, con tappo a vite blu, sterile R |
| 111547 | 15 ml Polypropylene test tube, Ø 17 x 120 mm, graduated up to 14 ml, with blue screw cap | Provetta in polipropilene da 15 ml, Ø 17 x 120 mm, graduata fino a 14 ml, con tappo a vite blu |
| 111548 | 15 ml Polypropylene test tube, Ø 17 x 120 mm, graduated up to 14 ml, with blue screw cap, sterile R | Provetta in polipropilene da 15 ml, Ø 17 x 120 mm, graduata fino a 14 ml, con tappo a vite blu, sterile R |
| 118548 | 15 ml Polypropylene test tube Ø 17 x 120 mm, graduated up to 14 ml, with blue screw cap, sterile, on expanded PS tray, 50 pcs, in sealed PE bag | Provetta in polipropilene da 15 ml, Ø 17 x 120 mm, graduata fino a 14 ml, con tappo a vite blu, su vassoio di Polistirolo espanso, in sacchetto di PE sigillato, sterile R |
| 119548 | 15 ml Polypropylene test tube Ø 17 x 120 mm, graduated up to 14 ml, with blue screw cap, sterile R individually wrapped | Provetta in polipropilene da 15 ml, Ø 17 x 120 mm, graduata fino a 14 ml, con tappo a vite blu, confezione singola, sterile R |
| 112013 | 12 ml Polystyrene conical test tube, Ø 16 x 110 mm, for Autochemyst AGA type | Provetta conica in polistirolo da 12 ml, Ø 16 x 110 mm, per autoanalizzatori tipo Auto Chemist AGA |
| 112017 | 12 ml Polystyrene urine test tube Ø 17 x 105 mm, graduated at 3-6-9-12 ml (to plug with blue stopper 113136) | Provetta in polistirolo da 12 ml, Ø 17 x 105 mm, per urina, graduata a 3-6-9-12 ml (da chiudere con tappo a pressione, blu, cod. 113136) |
| 113136 | Blue Polyethylene stopper Ø 17 mm for 112017 | Tappo blu in polietilene Ø 17 mm, per provetta urine 112017 |
| 112410 | Uriplast® 10 ml Polystyrene urine test tube with sediment collection receptacle, graduated at 2,5-5-10 ml (to plug with stopper 114137) | Provetta in polistirolo da 10 ml, Ø 16 x 105 mm, per urina Uriplast® con pozzetto di raccolta del sedimento, graduata a 2,5-5-10 ml (da chiudere con tappo a pressione 114137) |
| 114137 | Special hygienic natural Polyethylene stopper Ø 16 mm, for 112410 urine test tube | Tappo igienico in polietilene naturale Ø 16 mm, per provetta urina 112410 |
| 118710 | 8 ml Polystyrene cylindrical test tube Ø 13 x 100 mm, (to plug with yellow screw cap 118733) | Provetta cilindrica in polistirolo da 8 ml, Ø 13 x 100 mm, da chiudere con tappo a vite giallo 118733 |
| 115710 | 5 ml Polystyrene cylindrical test tube Ø 13 x 75 mm (to plug with yellow screw cap 118733) | Provetta cilindrica in polistirolo da 5 ml, Ø 13 x 75 mm, da chiudere con tappo a vite giallo 118733 |
| 118733 | Yellow Polyethylene stopper Ø 13 mm for 115710 and 118710 | Tappo giallo Ø 13 mm, in polietilene per provette 115710 e 118710 |
| 111047 | 10 ml Polypropylene gamma radiation resistant, cylindrical test tube, Ø 16 x 100 mm, graduated at 2,5 - 5 - 10 ml | Provetta cilindrica da 10 ml in polipropilene resistente alla sterilizzazione a raggi gamma, Ø 16 x 100 mm, graduata a 2,5 - 5 - 10 ml |
| 111040 | 10 ml Polypropylene cylindrical test tube, Ø 16 x 100 mm, graduated at 2,5 - 5 - 10 ml | Provetta cilindrica in polipropilene da 10 ml, Ø 16 x 100 mm, graduata a 2,5 - 5 - 10 ml |
| 111140 | 10 ml Polypropylene cylindrical test tube, Ø 16 x 100 mm, graduated at 2,5 - 5 - 10 ml, labelled | Provetta cilindrica in polipropilene da 10 ml, Ø 16 x 100 mm, graduata a 2,5 - 5 - 10 ml, etichettata |
| 111148 | 10 ml Polypropylene cylindrical test tube, Ø 16 x 100 mm, graduated at 2,5 - 5 - 10 ml, with yellow stopper, sterile R | Provetta cilindrica in polipropilene da 10 ml, Ø 16 x 100 mm, graduata a 2,5 - 5 - 10 ml, etichettata, tappo giallo, sterile R |
| 111010 | 10 ml Polystyrene cylindrical test tube, Ø 16 x 100 mm, graduated at 2,5 - 5 - 10 ml | Provetta cilindrica in polistirolo da 10 ml, Ø 16 x 100 mm, graduata a 2,5 - 5 - 10 ml |
| 111110 | 10 ml Polystyrene cylindrical test tube, Ø 16 x 100 mm, graduated at 2,5 - 5 - 10 ml, labelled | Provetta cilindrica in polistirolo da 10 ml, Ø 16 x 100 mm, graduata a 2,5 - 5 - 10 ml, etichettata singolarmente |
| 111018 | 10 ml Polystyrene cylindrical test tube, Ø 16 x 100 mm, graduated at 2,5 - 5 - 10 ml, sterile, plugged with yellow stopper | Provetta cilindrica in polistirolo da 10 ml, Ø 16 x 100 mm, graduata a 2,5 - 5 - 10 ml, con tappo giallo, sterile R |
| 111118 | 10 ml Clear Polystyrene cylindrical test tube, Ø 16 x 100 mm, graduated at 2,5 - 5 - 10 ml, sterile, plugged with yellow stopper, individually labelled | Provetta cilindrica da 10 ml in polistirolo, Ø 16 x 100 mm, graduata a 2,5 - 5 - 10 ml, con tappo giallo, etichettata singolarmente, sterile R |
| 112040 | 10 ml Polypropylene test tube, Ø 16 x 100 mm, graduated at 0,5- 1 - 2,5 - 5 - 10 ml | Provetta conica da 10 ml in polipropilene, Ø 16 x 100 mm, graduata a 0,5 - 1 - 2,5 - 5 - 10 ml |
| 112140 | 10 ml Polypropylene test tube, Ø 16 x 100 mm, graduated at 0,5-1- 2,5 - 5 - 10 ml, individually labelled | Provetta conica da 10 ml in polipropilene, Ø 16 x 100 mm, graduata a 0,5 - 1 - 2,5 - 5 - 10 ml, etichettata |





| Code | Description  | Descrizione  |
|--------|--|--|
| 112010 | 10 ml Polystyrene conical test tube, Ø 16 x 100 mm, graduated at 0,5 - 1 - 2,5 - 5 - 10 ml | Provetta conica da 10 ml in polistirolo, Ø 16 x 100 mm, graduata a 0,5 - 1 - 2,5 - 5 - 10 ml |
| 112110 | 10 ml Polystyrene conical test tube, Ø 16 x 100 mm, graduated at 0,5 - 1 - 2,5 - 5 - 10 ml, labelled | Provetta conica da 10 ml in polistirolo, Ø 16 x 100 mm, graduata a 0,5 - 1 - 2,5 - 5 - 10 ml, etichettata |
| 112018 | 10 ml Polystyrene conical test tube, Ø 16 x 100 mm, graduated at 0,5 - 1 - 2,5 - 5 - 10 ml, plugged with yellow stopper, sterile R | Provetta conica da 10 ml in polistirolo, Ø 16 x 100 mm, graduata a 0,5 - 1 - 2,5 - 5 - 10 ml, chiusa con tappo giallo, sterile R |
| 112118 | 10 ml Polystyrene conical test tube, Ø 16 x 100 mm, graduated at 0,5 - 1 - 2,5 - 5 - 10 ml, plugged with yellow stopper, labelled, sterile R | Provetta conica da 10 ml in polistirolo, Ø 16 x 100 mm, graduata a 0,5 - 1 - 2,5 - 5 - 10 ml, chiusa con tappo giallo, etichettata, sterile R |
| 112918 | 10 ml Polystyrene conical test tube, Ø 16 x 100 mm, graduated at 0,5 - 1 - 2,5 - 5 - 10 ml, plugged with yellow stopper, individually wrapped, sterile R | Provetta conica da 10 ml in Polistirolo trasparente, Ø 16 x 100 mm, graduata a 0,5 - 1 - 2,5 - 5 - 10 ml, chiusa con tappo giallo, confezione singola, sterile R |
| 119040 | 8 ml Polypropylene cylindrical test tube, Ø 16 x 75 mm, graduated at 2,5 - 5 ml | Provetta cilindrica da 8 ml in polipropilene, Ø 16 x 75 mm graduata a 2,5 - 5 ml |
| 119140 | 8 ml Polypropylene cylindrical test tube, Ø 16 x 75 mm, graduated at 2,5 - 5 ml, labelled | Provetta cilindrica da 8 ml in polipropilene, Ø 16 x 75 mm, graduata a 2,5 - 5 ml, etichettata |
| 119010 | 8 ml Polystyrene cylindrical test tube, Ø 16 x 75 mm, graduated at 2,5 - 5 ml | Provetta cilindrica da 8 ml in polistirolo, Ø 16 x 75 mm, graduata a 2,5 - 5 ml |
| 119110 | 8 ml Polystyrene cylindrical test tube, Ø 16 x 75 mm, graduated at 2,5 - 5 ml, individually labelled | Provetta cilindrica da 8 ml in polistirolo, Ø 16 x 75 mm, graduata a 2,5 - 5 ml, etichettata |
| 113040 | 5 ml Polypropylene cylindrical test tube, Ø 12 x 92mm graduated at 0,5 - 2,5 - 5 ml | Provetta cilindrica da 5 ml in polipropilene, Ø 12 x 92 mm, graduata a 0,5 - 2,5 - 5 ml |
| 113140 | 5 ml Polypropylene cylindrical test tube, Ø 12 x 92mm graduated at 0,5 - 2,5 - 5 ml, individually labelled | Provetta cilindrica da 5 ml in polipropilene, Ø 12 x 92 mm, graduata a 0,5 - 2,5 - 5 ml, etichettata |
| 113010 | 5 ml Polystyrene cylindrical test tube, Ø 12 x 92 mm graduated at 0,5 - 2,5 - 5 ml | Provetta cilindrica da 5 ml, in polistirolo, Ø 12 x 92 mm graduata a 0,5 - 2,5 - 5 ml |
| 113110 | 5 ml Polystyrene cylindrical test tube, Ø 12 x 92 mm graduated at 0,5 - 2,5 - 5 ml, labelled | Provetta cilindrica da 5 ml, in polistirolo, Ø 12 x 92 mm, graduata a 0,5 - 2,5 - 5 ml, etichettata |
| 113018 | 5 ml Polystyrene cylindrical test tube, Ø 12 x 92 mm graduated at 0,5 - 2,5 - 5 ml, plugged with yellow stopper, sterile R | Provetta cilindrica da 5 ml, in polistirolo, Ø 12 x 92 mm, graduata a 0,5 - 2,5 - 5 ml, chiusa con tappo giallo, sterile R |
| 113118 | 5 ml Polystyrene cylindrical test tube, Ø 12 x 92 mm graduated at 0,5 - 2,5 - 5 ml, plugged with yellow stopper, labelled, sterile R | Provetta cilindrica da 5 ml, in polistirolo, Ø 12 x 92 mm, graduata a 0,5 - 2,5 - 5 ml, chiusa con tappo giallo, etichettata, sterile R |
| 110440 | 5 ml Polypropylene cylindrical test tube, Ø 12 x 75 mm, with rim | Provetta cilindrica da 5 ml, in polipropilene Ø 12 x 75 mm, con bordo |
| 110441 | 5 ml Polypropylene cylindrical test tube, Ø 12 x 75 mm, without rim, | Provetta cilindrica da 5 ml in polipropilene senza bordo, Ø 12 x 75 mm |
| 110442 | 5 ml Polypropylene cylindrical test tube, Ø 12 x 75 mm, without rim, 250 pcs per bag | Provetta cilindrica da 5 ml in polipropilene Ø 12 x 75 mm, senza bordo, in sacchetti da 250 pz. |
| 110411 | 5 ml Polystyrene cylindrical test tube without rim, Ø 12 x 75 mm | Provetta cilindrica da 5 ml in polistirolo, Ø 12 x 75 mm, senza bordo |
| 110412 | 5 ml Clear Polystyrene cylindrical test tube without rim, Ø 12 x 75 mm - 250 pcs per bag | Provetta cilindrica da 5 ml in polistirolo, Ø 12 x 75 mm, senza bordo, in sacchetti da 250 pz. |
| 110418 | 5 ml Polystyrene cylindrical test tube, Ø 12 x 75 mm, without rim, with yellow stopper, in bags of 250 pcs, sterile R | Provetta cilindrica da 5 ml in polistirolo, Ø 12 x 75 mm, senza bordo, chiusa con tappo giallo, in sacchetti da 250 pz, sterile R |
| 115040 | 5 ml Polypropylene cylindrical test tube Ø 16 x 60mm, flat bottom, graduated at 1 - 2,5 - 5 ml | Provetta cilindrica da 5 ml in polipropilene, Ø 16 x 60 mm, fondo piatto, graduata a 1 - 2,5 - 5 ml |
| 115140 | 5 ml Polypropylene cylindrical test tube Ø 16 x 60mm, flat bottom, graduated at 1 - 2,5 - 5 ml, labelled | Provetta cilindrica da 5 ml in polipropilene, Ø 16 x 60 mm, fondo piatto, graduata a 1 - 2,5 - 5 ml, etichettata |
| 115010 | 5 ml Polystyrene cylindrical test tube Ø 16 x 60 mm, flat bottom, graduated at 1 - 2,5 - 5 ml | Provetta cilindrica da 5 ml in polistirolo, Ø 16 x 60 mm, fondo piatto, graduata a 1 - 2,5 - 5 ml |
| 115110 | 5 ml Polystyrene cylindrical test tube Ø 16 x 60 mm, flat bottom, graduated at 1 - 2,5 - 5 ml, labelled | Provetta cilindrica da 5 ml in polistirolo, Ø 16 x 60 mm, fondo piatto, graduata a 1 - 2,5 - 5 ml, etichettata |
| 115048 | 5 ml Polypropylene cylindrical test tube Ø 16 x 60mm, flat bottom, graduated at 1 - 2,5 - 5 ml sterile, plugged with yellow stopper | Provetta cilindrica da 5 ml in polipropilene, fondo piatto, Ø 16 x 60 mm, graduata a 1 - 2,5 - 5 ml, chiusa con tappo giallo, sterile R |
| 115148 | 5 ml Polypropylene cylindrical test tube Ø 16 x 60mm, flat bottom, graduated at 1 - 2,5 - 5 ml sterile, plugged with yellow stopper, individually labelled | Provetta cilindrica da 5 ml in polipropilene, fondo piatto, Ø 16 x 60 mm, graduata a 1 - 2,5 - 5 ml, chiusa con tappo giallo, etichettata, sterile R |
| 110510 | 4 ml Polystyrene cylindrical test tube Ø 13 x 50 mm, graduated at 1 - 2,5 - 5 ml | Provetta cilindrica da 4 ml in polistirolo, Ø 13 x 50 mm, graduata a 1 - 2,5 - 5 ml |



| Code | Description  | Descrizione  |
|--------|--|--|
| 117041 | 3 ml Polypropylene cylindrical test tube w/out rim, Ø 11,5 x 55 mm | Provetta cilindrica da 3 ml in polipropilene, Ø 11,5 x 55 mm, senza bordo |
| 117011 | 3 ml Polystyrene cylindrical test tube w/out rim, Ø 11,5 x 55 mm | Provetta cilindrica da 3 ml in polistirolo, Ø 11,5 x 55 mm, senza bordo |
| 118011 | 3 ml Polystyrene cylindrical test tube Ø 11 x 70 mm for RIA, graduated at 2 ml, to plug with stopper 112430 | Provetta cilindrica da 3 ml in polistirolo, Ø 11 x 70 mm, in polistirolo, per RIA, graduata a 2 ml, da chiudere con tappo cod. 112430 |
| 112430 | Natural Polyethylene stopper Ø 11 mm, for 118011 test tube | Tappo in polietilene naturale a presa alta, Ø 11 mm, per provetta 118011 |
| 117040 | 3 ml Polypropylene cylindrical test tube Ø 12 x 56 mm graduated at 0,5 - 1 - 2 - 3 ml | Provetta cilindrica da 3 ml in polipropilene, Ø 12 x 56 mm, in polipropilene, graduata a 0,5 - 1 - 2 - 3 ml |
| 117140 | 3 ml Polypropylene cylindrical test tube Ø 12 x 56 mm graduated at 0,5 - 1 - 2 - 3 ml, individually labelled | Provetta cilindrica da 3 ml in polipropilene, Ø 12 x 56 mm, in polipropilene, graduata a 0,5 - 1 - 2 - 3 ml, etichettata singolarmente |
| 117010 | 3 ml Polystyrene cylindrical test tube Ø 12 x 56 mm, graduated at 0,5 - 1 - 2 - 3 ml | Provetta cilindrica da 3 ml in polistirolo, Ø 12 x 56 mm, in polistirolo, graduata a 0,5 - 1 - 2 - 3 ml |
| 117110 | 3 ml Polystyrene cylindrical test tube Ø 12 x 56 mm, graduated at 0,5 - 1 - 2 - 3 ml individually labelled | Provetta cilindrica da 3 ml in polistirolo, Ø 12 x 56 mm, in polistirolo, graduata a 0,5 - 1 - 2 - 3 ml, etichettata |
| 117018 | 3 ml Polystyrene cylindrical test tube Ø 12 x 56 mm, graduated at 0,5 - 1 - 2 - 3 ml sterile, plugged with yellow stopper | Provetta cilindrica da 3 ml in polistirolo, Ø 12 x 56 mm, in polistirolo, graduata a 0,5 - 1 - 2 - 3 ml, chiusa con tappo giallo, sterile R |
| 117118 | 3 ml Polystyrene cylindrical test tube Ø 12 x 56 mm, graduated at 0,5 - 1 - 2 - 3 ml sterile, plugged with yellow stopper labelled | Provetta cilindrica da 3 ml in polistirolo, Ø 12 x 56 mm, graduata a 0,5 - 1 - 2 - 3 ml, chiusa con tappo giallo, etichettata, sterile R |
| 114040 | 3 ml Polypropylene cylindrical test tube Ø 12 x 51mm, flat bottom, graduated at 0,5 - 1 - 2 - 2,5 ml | Provetta cilindrica da 3 ml in polipropilene, Ø 12 x 51 mm, fondo piatto, graduata a 0,5 - 1 - 2 - 2,5 ml |
| 114140 | 3 ml Polystyrene cylindrical test tube Ø 12 x 51 mm, flat bottom, graduated at 0,5 - 1 - 2 - 2,5 ml, labelled | Provetta cilindrica da 3 ml in polipropilene, Ø 12 x 51 mm, fondo piatto, graduata a 0,5 - 1 - 2 - 2,5 ml, etichettata |
| 114010 | 3 ml Polystyrene cylindrical test tube Ø 12 x 51 mm, flat bottom, graduated at 0,5 - 1 - 2 - 2,5 ml | Provetta cilindrica da 3 ml in polistirolo, Ø 12 x 51 mm, fondo piatto, graduata a 0,5 - 1 - 2 - 2,5 ml |
| 114110 | 3 ml Polystyrene cylindrical test tube Ø 12 x 51 mm, flat bottom, graduated at 0,5 - 1 - 2 - 2,5 ml, labelled | Provetta cilindrica da 3 ml in polistirolo, Ø 12 x 51 mm, fondo piatto, graduata a 0,5 - 1 - 2 - 2,5 ml, etichettata |
| 114048 | 3 ml Polypropylene cylindrical test tube Ø 12 x 51mm, flat bottom, graduated at 0,5 - 1 - 2 - 2,5 ml, plugged with yellow stopper, sterile R | Provetta cilindrica da 3 ml in polipropilene, Ø 12 x 51 mm, fondo piatto, graduata a 0,5 - 1 - 2 - 2,5 ml, chiusa con tappo giallo, sterile R |
| 114148 | 3 ml Polypropylene cylindrical test tube Ø 12 x 51mm, flat bottom, graduated at 0,5 - 1 - 2 - 2,5 ml, plugged with yellow stopper, labelled, sterile R | Provetta cilindrica da 3 ml in polipropilene, Ø 12 x 51 mm, fondo piatto, graduata a 0,5 - 1 - 2 - 2,5 ml, chiusa con tappo giallo, etichettata, sterile R |
| 110041 | 1 ml Polypropylene cylindrical test tube Ø 8 x 40 mm, without rim | Provetta cilindrica da 1 ml, in polipropilene, Ø 8 x 40 mm, senza bordo |
| 110110 | 20 ml Polystyrene cylindrical test tube without rim, Ø 16 x 150 mm | Provetta cilindrica da 20 ml in polistirolo, Ø 16 x 150 mm, senza bordo, |
| 110118 | 20 ml Polystyrene cylindrical test tube without rim, Ø 16 x 150 mm sterile, with yellow stopper | Provetta cilindrica da 20 ml in polistirolo, Ø 16 x 150 mm, senza bordo, con tappo giallo, sterile R |
| 110210 | 15 ml Polystyrene cylindrical test tube without rim, Ø 16 x 125 mm | Provetta cilindrica da 15 ml in polistirolo, Ø 16 x 125 mm, senza bordo, |
| 261170 | Venoplast® tube Ø 10 x 75 mm for ESR by Sediten® for 1.6 ml of blood, no anticoagulant, no vacuum, in peel pack of 100 pcs. | Provetta Venoplast® Ø 10 x 75 mm per VES su Sediten®, per 1,6 ml di sangue, senza vuoto, senza anticoagulante, in blister da 100 pz |
| 290300 | 7 ml cylindrical test tube Ø 12 x 100 mm, in transparent TPX, 170°C heat resistant, autoclavable | Provetta cilindrica da 7 ml in TPX, Ø 12 x 100 mm, resistente al calore fino a 170°C, autoclavabile |
| 290301 | 10 ml cylindrical test tube Ø 16 x 100 mm, in transparent TPX, 170°C heat resistant, autoclavable | Provetta cilindrica da 10 ml in TPX, Ø 16 x 100 mm, resistente al calore fino a 170°C, autoclavabile |
| 290302 | 16 ml cylindrical test tube Ø 18 x 100 mm, in transparent TPX, 170°C heat resistant, autoclavable | Provetta cilindrica da 16 ml in TPX, Ø 18 x 100 mm, resistente al calore fino a 170°C, autoclavabile |
| 290303 | 26 ml cylindrical test tube Ø 24 x 90 mm, in transparent TPX, 170°C heat resistant, autoclavable | Provetta cilindrica da 26 ml in TPX, Ø 24 x 90 mm, resistente al calore fino a 170°C, autoclavabile |
| 290304 | 30 ml cylindrical test tube, Ø 20 x 100 mm, in transparent TPX, 170°C heat resistant, autoclavable | Provetta cilindrica da 30 ml in TPX, Ø 20 x 100 mm, resistente al calore fino a 170°C, autoclavabile |
| 290305 | 31 ml cylindrical test tube, Ø 24 x 95 mm, in Transparent TPX, heat resistant, autoclavable | Provetta cilindrica da 31 ml, in TPX, Ø 24 x 95 mm, resistente al calore fino a 170°C, autoclavabile |
| 290306 | 48 ml cylindrical test tube, Ø 30 x 100 mm, in Transparent TPX, heat resistant, autoclavable | Provetta cilindrica da 48 ml, in TPX, Ø 30 x 100 mm, resistente al calore fino a 170°C, autoclavabile |
| 290307 | 75 ml cylindrical test tube, Ø 35 x 100 mm, in Transparent TPX, heat resistant, autoclavable | Provetta cilindrica da 75 ml, in TPX Ø 35 x 100 mm, resistente al calore fino a 170°C, autoclavabile |
| 290308 | 110 ml cylindrical test tube, Ø 40 x 120 mm, in Transparent TPX, heat resistant, autoclavable | Provetta cilindrica da 110 ml, in TPX Ø 40 x 120 mm, resistente al calore fino a 170°C, autoclavabile |



| Code | Description  | Descrizione  |
|--------|--|---|
| 290309 | 160 ml cylindrical test tube, Ø 45 x 120 mm, in Transparent TPX, heat resistant, autoclavable | Provetta cilindrica da 160 ml, in TPX Ø 45 x 120 mm, resistente al calore fino a 170°C, autoclavabile |
| 290310 | 200 ml cylindrical test tube, Ø 50 x 120 mm, in Transparent TPX, heat resistant, autoclavable | Provetta cilindrica da 200 ml, in TPX Ø 50 x 120 mm, resistente al calore fino a 170°C, autoclavabile |
| 290320 | 7 ml cylindrical test tube, Ø 12 x 100 mm, in autoclavable Polypropylene, high chemical resistance | Provetta cilindrica da 7 ml, Ø 12 x 100 mm, in polipropilene autoclavabile a 121°C resistente agli acidi |
| 290321 | 10 ml cylindrical test tube, Ø 16 x 100 mm, in autoclavable Polypropylene, high chemical resistance | Provetta cilindrica da 10 ml, Ø 16 x 100 mm, in polipropilene autoclavabile a 121°C resistente agli acidi |
| 290322 | 16 ml cylindrical test tube, Ø 18 x 100 mm, in autoclavable Polypropylene, high chemical resistance | Provetta cilindrica da 16 ml, Ø 18 x 100 mm, in polipropilene autoclavabile a 121°C resistente agli acidi |
| 290323 | 26 ml cylindrical test tube, Ø 24 x 90 mm, in autoclavable Polypropylene, high chemical resistance | Provetta cilindrica da 26 ml, Ø 24 x 90 mm, in polipropilene autoclavabile a 121°C resistente agli acidi |
| 290324 | 30 ml cylindrical test tube, Ø 20 x 100 mm, in autoclavable Polypropylene, high chemical resistance | Provetta cilindrica da 30 ml, Ø 20 x 100 mm, in polipropilene autoclavabile a 121°C, resistente agli acidi |
| 290325 | 31 ml cylindrical test tube, Ø 24 x 95 mm, in autoclavable Polypropylene, high chemical resistance | Provetta cilindrica da 31 ml, Ø 24 x 95 mm, in polipropilene autoclavabile a 121°C, resistente agli acidi |
| 290326 | 48 ml cylindrical test tube, Ø 30 x 100 mm, in autoclavable Polypropylene, high chemical resistance | Provetta cilindrica da 48 ml, Ø 30 x 100 mm, in polipropilene autoclavabile a 121°C, resistente agli acidi |
| 290327 | 75 ml cylindrical test tube, Ø 35 x 100 mm, in autoclavable Polypropylene, high chemical resistance | Provetta cilindrica da 75 ml, Ø 35 x 100 mm, in polipropilene autoclavabile a 121°C, resistente agli acidi |
| 290328 | 110 ml cylindrical test tube, Ø 40 x 120 mm, in autoclavable Polypropylene, high chemical resistance | Provetta cilindrica da 110 ml, Ø 40 x 120 mm, in polipropilene autoclavabile a 121°C, resistente agli acidi |
| 290329 | 160 ml cylindrical test tube, Ø 45 x 120 mm, in autoclavable Polypropylene, high chemical resistance | Provetta cilindrica da 160 ml, Ø 45 x 120 mm, in polipropilene autoclavabile a 121°C, resistente agli acidi |
| 290330 | 200 ml cylindrical test tube, Ø 50 x 120 mm, in autoclavable Polypropylene, high chemical resistance | Provetta cilindrica da 200 ml, Ø 50 x 120 mm, in polipropilene autoclavabile a 121°C, resistente agli acidi |
| 290335 | 10 ml conical test tube, Ø 16 x 108 mm, in Transparent TPX, heat resistant, autoclavable | Provetta conica da 10 ml, Ø 16 x 108 mm, in TPX, autoclavabile, resistente al calore fino a 170°C |
| 290336 | 16 ml conical test tube, Ø 18 x 120 mm, in Transparent TPX, heat resistant, autoclavable | Provetta conica da 16 ml, Ø 18 x 120 mm, in TPX autoclavabile, resistente al calore fino a 170°C |
| 290340 | 10 ml conical test tube, Ø 16 x 108 mm, in autoclavable Polypropylene, high chemical resistance | Provetta conica da 10 ml, Ø 16 x 108 mm, in polipropilene autoclavabile a 121°C, resistente agli acidi |
| 290341 | 16 ml conical test tube, Ø 18 x 120 mm, in autoclavable Polypropylene, high chemical resistance | Provetta conica da 16 ml, Ø 18 x 120 mm, in polipropilene autoclavabile a 121°C, resistente agli acidi |